


Interlig + Fibrex-Lab Pontic



With these products you can guarantee immediate esthetics and strength.

Interlig and the Pontic System allow the fabrication of immediate resin-retained fixed bridges which are highly resistant, esthetic and long-lasting.



1. Initial case.

2. Select the Pontic size that best fits the prosthetic space.

3. Insert Interlig on the groove and on dental preparations, adapt it with the fluid resin.

4. Light cure the structure.

5. Build the bridge with composite resin.

Clinical case:
Prof. Raphael Monte Alto



CHANGE LANGUAGE